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### 1. Instructions for use and safety, important information

All products of Denti Implant Systems are recommended exclusively for dentists and oral surgeons who are skilled in dental implantology and oral surgery. In every case, it is the responsibility of the user to decide about the proper utilization of Denti Implant Systems in accordance with proper medical examinations. It is also the responsibility of the user of Denti Implant Systems to decide whether the products are appropriate for a certain patient in a given case. Denti System Ltd. disclaims any obligations or responsibility for any direct or indirect damage caused by the implantation or utilization of Denti products among conditions that resulted from decisions which were not appropriately made or which were the results of inexperience of the user.

All elements of this surgical and prosthesis set are part of an integrated system, therefore this set can only be utilized safely and successfully – in accordance with instructions and guides included by the supplier – with its original components and devices. In any case any part of these sets is used along with any foreign device, this may negatively influence the results, which exempts Denti System Ltd. from any guarantee or other events involving legal procedures. Every device and tool within Denti Implant Systems has been developed for special applications, any incorrect utilization may cause these devices to wear off sooner than usual and thus endanger the success of application.

Drills, devices and all other components of Denti implant Systems have been developed for certain types and diameters. Utilization with inappropriate types or with inappropriate diameters may cause mechanical damage to the components of the system or even tissue damage to the patient. There is an appropriate surgical set for each implant type. Color marks help with the selection of necessary devices and components.

Studying the details of the utilization of Denti Implant Systems is the responsibility of the user. In case of any doubt or dubious situations, the user is obliged to contact Denti System Ltd.

### I. IMPLANTS

#### 2. Features and details of Denti implants

Denti implants are made of high purity titanium designed for use in humans (Grade 2 and Grade 4 quality titanium). The surface of the implants is first sandblasted with large granule bioceramics, and then etched with acid - or in other words, passivated. Abutments and assemblies are made of Ti6Al14V titanium compounds, circonium oxide ceramics and polymer materials.

#### 3. Indications

Denti Implant System implants can be utilized to provide implantation prosthesis for the functional or esthetic rehabilitation of both the lower and the upper jaws. Implants can be used in case of short defects in the middle

or the row or even to compensate for partial or even complete absences of denture. Head parts can be used to make prostheses fixed by means of cementing or screwing or even to make removable prostheses. All indications of various types of Denti implants and head parts as well as instructions for special use can be found in guides supplied with Denti Implant Systems.

Additional indications

*Denti Cylindric Screw, Root Form and Bone-Level Implants* can be implanted during surgical exploration (sub-gingivally - thus closed when healing) or transgingivally (thus being open when healing), or during a one- or two-phase surgery when implanted into both the upper and the lower jaws.

*Denti OP and Needle one-stage root-shaped tapered implants* can be implanted into both the upper and the lower jaws by means of the single-stage surgical method. Denti OP and Denti N implants can be used for both immediate and delayed loading. In case bone structures are sufficient both in quantity and in quality, fixed implantation prostheses, which have been designed and made for the immediate or early functional loading of OP and Denti N implants, can be made upon patient request. These implants can be used as stand-alone elements to fix supported crowns and bridges. They can also be used to make so-called mixedly supported prostheses with teeth having a resilient periodontium. They can also be applied to make temporary implants between traditionally loaded implants so that a temporary prosthesis can be worn during the period of osseointegration.

These implants can be loaded for any indications (as oral implants to be loaded into the bone (a single missing tooth, a bridge, partial or complete prosthesis)) or as extraoral implants (see Denti Extraoral implants) to compensate for any defects as part of restoring or esthetic surgery (eye, nose, ear epitheses and protheses). In case the bone tissue is sufficient, we recommend using the cylindric screw implant with a diameter of 4. 0-4. 5 Ø and root-shaped implants with a diameter of 4. 3-4. 8 Ø. For any further details of implantation see Denti manuals.

#### 4. Contraindications

The most important contraindications of implantation are the following: recurring mucosal conditions, inflammations of the jaws and gingiva, poor or insufficient oral hygiene, dryness of the mouth (xerostomia), untreated parafunctional habits (scrunching, lockjaw), smoking, and being under 18 years of age. Implantations are also contraindicated in case of any alterations or conditions that may complicate surgery, including neurological and psychological diseases, disorders or any other condition that may hinder the patient's practice of regular oral hygiene. Additionally, implants are contraindicated in case of any bone alterations due to metabolic disorders, untreated hemophilia, or tumors, and after the treatment of tumors in the next 1. 5-2 years. Other contraindicating systemic diseases include blood dyscrasia, untreated diabetes, hyperthyroidism, AIDS, drug abuse, excessive alcohol consumption or osteoporosis. Implantations are also contraindicated in case the patient had a myocardial infarction in the previous 12 months, or in case the patient is immunocompromised or has leukocytosis, or taking steroids regularly, having titanium allergy, or untreated endocrine disorders.

<div><div><span><span> </span></span></div><span> </span></div>	3,3 and 3. 4 <span> </span> mm implants should not be loaded to make up for the single upper first incisors, or any single canine, premolar or molar teeth.
---	---

Loading OP or Denti N implants into soft (D-4) bone tissue is not recommended. In case of D-3 bone quality, please consider implantation carefully. In case there is a circumscribed, significant bone deficiency at the site of implantation, or in case the bone is not sufficient to take in the entire implant, then implantation is again not recommended. Loading Denti OP implants is also not recommended in case the vertical distance between the antagonist chewing surface and the alveolar crest is less than 8 mm. Lateral and/or vertical augmentation should not be more than 2 mm for OP and N implants. Denti OP and Denti N conical head-part implants can only support removable prosthesis as telescopic fixation.

#### 5. Side effects

After surgical implantation, patients should refrain from any physical stress and smoking. Medical instructions and the prescribed therapy must be followed carefully!
Temporary complaints may include pain, facial swelling, speech disorders, gingivitis, hematomas or late onset bleeding. Long-term complaints may include permanent pain caused by the dental implant, permanent par- esthesia, dysesthesia, loss of bone tissue from the maxillary/mandibular ridge, local or systemic infections / oroantral or oronasal inflammations/, loss of neighboring teeth, irreversible damage to the neighboring teeth, esthetic problems, local or systemic allergic response. It is the responsibility of the dentist to select patients who are expected to be satisfied psychologically, functionally and esthetically after recovery.

#### 6. Warning

Implantation surgery and implantation prosthetics are complex dental procedures. It is especially recommended to master all necessary elements of proficiency. An improper surgical procedure may result in implant loss or bone degeneration. The implantation must be preceded by extensive diagnostic procedures. Improper design may result in unsuccessful implantation. Denti implants must only be loaded with devices included in the Denti Implant Systems and in accordance with guidelines and recommendations supplied therewith. Studying the instructions for use alone is not enough to start utilizing Denti Implant Systems. Denti implants may exclusively be loaded by oral surgeons, dentists and specialists trained to use Denti Systems. Denti System Ltd. holds regular trainings, which allow for a deeper knowledge in this field. Guarantees shall be submitted only if the claimer is able to confirm the completion of a Denti System Ltd. course.

#### 7. Sterilization and packing

Implants are released following gamma sterilization and are made for a single use only. Clinically infected implants must not be cleaned and resterilized. Denti System Ltd. does not take any liability for the utilization of Denti implants taken from open, resterilized Denti implant sets. Used or unsterile implants must not be loaded. The sterile packing protects Denti implants from outer contamination. Inner packing must only be opened directly before loading. Prior to loading, personnel should check that the packing was not damaged in a way that the implant sterility was compromised. When stored properly (at a dry place, protected from sunlight at room temperature and in its original packing) implants remain sterile till the date of expiry printed on the packing. Implants must not be used after the date of expiry!

All Denti implants can be traced back according to product name and registration no., which can be found on both the outer and the inner packings.

#### 8. Handling the implants

Denti *titanium* implants are packed in two-fold packing, thus the surface of implants is protected from damage or contamination. After removing the protective packing, it is strictly forbidden to touch the implants by hand

or any other way. The osseointegration surface of the implants must be touched with tools made from pure titanium only!

#### PREPARING THE IMPLANTS

Each and every Denti implant is packed in a double-layer sterile packing. Before use, implants should be stored at room temperature, at a dry place, protected from sunlight. Before opening, check the date of expiry and whether there is any damage on the packing. Packing should be opened directly before use only. Implants must not be used if the date expired, or the packing was damaged before use or the packing is opened.

Elements and tools for the implants are quite small, so pay special attention to prevent the patient from accidentally swallowing or aspiring them.

Please mind all recommendations in Product Catalogues issued by Denti System Ltd. and recommendations found in guides and publications of Denti Implant Systems.

#### 9. Exchange

<div><div><span><span> </span></span></div><span> </span></div>	In case the original packing is damaged, Denti System Ltd. cannot exchange any implants.
---	--

#### 10. Validity

All previous versions of this guide are hereby made invalid by this version.

#### 11. Packing units

In case there are no other specifications involved, implants are packed one by one.

Failure to take any of the instructions detailed above excludes any and all liability of the company.

#### 12. Signs and explanations

The signs for tools fitted with color and numerical codes can be found on the plastic handle.

**Color codes**

Yellow for 3.3-3. 4-3. 8 Ø implants
Red for 4. 0-4. 3-4. 2. 3 Ø implants
Blue for 4. 5-4. 8-2. 8 Ø implants
Green for 5. 5-5. 3 Ø Ø implants
Black general tools
White 4. 0-4. 3-4. 5-4. 8-5. 3-5. 5 Ø implants

In order to identify the correct tools and sizes please use Denti Implant System Product Catalogue.

#### 13. Documentations and guarantee

Denti System Ltd recommends its customers to preserve the entire documentation of single cases. The inner packing of products contains the product code and registry number. Inner packing also contains a label with the product number, LOT (individual registry) number and a bar code of the product which is recommended to be documented along with patient medical data, so that implants can be identified when necessary. The guarantee of any product is only valid when LOT numbers and bar codes are both available.

For technical details or further information please refer to your local representative.

#### 14. Notation

Registry number – individual product identification

CE markings

Date of manufacture (year, month)

Do not use more than once

Product code

This product was sterilized with gamma rays

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