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INSTRUCTIONS FOR USAGE AND SAFETY

Instructions for Denti® Implant Systems

Denti® Root Form Implants (DR)
Denti® Root Form+ Implants (DR+)
Denti® One Piece Implants (DOP)
Denti® Bone Level Implants (DBL)
Denti® Kone Implants (DK)
Denti® Pro Implants (DP)
Denti® Needle Form Implants (DN, DNB)
Denti® Stabile 1 Implants (DF)

1. Instructions for usage and safety, important information

All products of Denti Implant Systems are recommended exclusively for dentists and oral surgeons who are skilled in dental implantology and oral surgery. In every case, it is the responsibility of the user to decide about the proper utilization of Denti Implant Systems in accordance with proper medical examinations. It is also the responsibility of the user of Denti Implant Systems to decide whether the products are appropriate for a certain patient in a given case. Denti System Ltd. disclaims any obligations or responsibility for any direct or indirect damage caused by the implantation or utilization of Denti products among conditions that resulted from decisions which were not appropriately made or which were the results of inexperience of the user. All elements of this surgical and prosthesis set are part of an integrated system, therefore this set can only be utilized safely and successfully – in accordance with instructions and guides included by the supplier – with its original components and devices. In case any part of these sets is used along with any foreign device, this may negatively influence the results, which exempts Denti System Ltd. from any guarantee or other events involving legal procedures. Every device and tool within Denti Implant Systems has been developed for special applications, any incorrect utilization may cause these devices to wear off sooner than usual and thus endanger the success of application.

Drills, devices and all other components of Denti Implant Systems have been developed for certain types and diameters. Utilization with inappropriate types or with inappropriate diameters may cause mechanical damage to the components of the system or even tissue damage to the patient. There is an appropriate surgical set for each implant type. Color marks help with the selection of necessary devices and components.

Studying the details of the utilization of Denti Implant Systems is the responsibility of the user. In case of any doubt or dubious situations, the user is obliged to contact Denti System Ltd.

I. IMPLANTS

2. Features and details of Denti implants

Denti implants are made of high purity titanium designed for use in humans (Grade 4 quality titanium). The surface of the implants is first sandblasted with large granule bioceramics, and then etched with acid - in other words, passivated. Abutments and assemblies are made of (Grade 5 quality titanium) Ti6Al4V ELI titanium compounds, chromium oxide ceramics and polymer materials.

3. Indications

Denti Implant System implants can be utilized to provide implantation prosthesis for the functional or esthetic rehabilitation of both the lower and the upper jaws. Implants can be used in case of short defects in the middle of the row or even to compensate for partial or even complete absences of denture. Head parts can be used to make prostheses fixed by means of cementing or screwing or even to make removable prostheses. All indications of various types of Denti implants and head parts as well as instructions for special use can be found in guides supplied with Denti Implant Systems.

Additional indications

Denti DR, DR+, DBL, DK, DP can be implanted during surgical exploration (subgingivally - thus closed when healing) or transgingivally (thus being open when healing), or during a one- or two-phase surgery when implanted into both the upper and the lower jaws.

Denti DOP, DN, DNB, DF one-stage root-shaped tapered implants can be implanted into both the upper and the lower jaws by means of the single-stage surgical method. Denti DOP, DN, DNB, DF implants can be used for both immediate and delayed loading. In case bone structures are sufficient both in quantity and in quality, fixed implantation prostheses, which have been designed and made for the immediate or early functional loading of DOP, DN, DNB, DF implants, can be made upon patient request. They can also be used to make so-called mixedly supported prostheses with teeth having a resilient periodontium. They can also be applied to make temporary implants between traditionally loaded implants so that a temporary prosthesis can be worn during the period of osseointegration.

For any further details of implantation see Denti manuals.

4. Contraindications

The most important contraindications of implantation are the following: recurring mucosal conditions, inflammations of the jaws and gingiva, poor or insufficient oral hygiene, dryness of the mouth (xerostomia), untreated parafunctional habits (scrunching, grinding lockjaw), smoking, and still growing jaws (typically under the age of 18).

Dental implant treatment is also contraindicated in case of any alterations or conditions that may complicate surgery, including neurological and psychological diseases, disorders or any other condition that may hinder the patient's practice of regular oral hygiene. Additionally, dental implant treatment are contraindicated in case of any bone alterations due to metabolic disorders (bone metabolic disorders, any disaese influencing bone regeneration or micro-circulation, rheumatic disease), untreated hemophilia, or tumors, and after the treatment of tumors in the next 1. 5-2 years. Other contraindicating systemic diseases include blood dyscrasia, haemophilia, untreated diabetes, hyperthyroidism, AIDS, drug abuse, excessive alcohol consumption or osteoporosis, serious systematic diseases, decreased immunity that increasesinfection risk, leukocyte dysfunction, endocrine disease, anticoagulation treatment, hemorrhagic diathesis arteriosclerosis, apoplexy, hypertension, hepatitis, diabetes mellitus. Implantations are also contraindicated in case the patient had a myocardial infarction in the previous

12 months, or in case the patient is immunocompromised or has leukocytosis, or taking steroids regularly, including titanium alloy, or untreated endocrine disorders. Inappropriate bone quality, radix relicta in the jaw, bone disease or insufficient bone and soft tissue quantity, healing difficulties, oral infections, serious therapy resistant functional disorder, Immunsuppressive treatment, connective tissue diseases/collagenosis, malignoma. Unfinished bone growth, anatomically insufficient bone status, jaw bone radiation therapy, tempo-mandibular joint disease (TMJ dysfunction), untreated jaw diseases.

Dental implant treatment is contraindicated for patients who treated by bisphosphonate and/ or denosumab containing medicines.

Important information:

3,3 mm, 3,4 mm and 3,8 mm implants should not be loaded to make up for the single upper first incisors, or any single canine, premolar or molar teeth.

One-phase implants should not be used for replacement of an individual tooth!

Loading DOP, DN, DNB and DF implants into soft (D-4) bone tissue is not recommended. In case of D-3 bone quality, please consider implantation carefully. In case there is a circumscribed, significant bone deficiency at the site of implantation, or in case the bone is not sufficient to take in the entire implant, then implantation is again not recommended. Loading DOP implants is also not recommended in case the vertical distance between the antagonist chewing surface and the alveolar crest is less than 8 mm. Lateral and/or vertical augmentation should not be more than 2 mm for DOP and NI implants. DOP and conical head-part implants can only support removable prosthesis as telescopic fixation.

5. Side effects

After surgical implantation, patients should refrain from any physical stress and smoking. Medical instructions and the prescribed therapy must be followed carefully!

Temporary complaints may include pain, facial swelling, speech disorders, gingivitis, hematomas or late onset bleeding. Long-term complaints may include permanent pain caused by the dental implant, permanent par-esthesia, dysesthesia, plate peeling, loss of bone tissue from the maxillary/mandibular ridge, local or systemic infections /ororantal or oronasal inflammations/, loss of neighboring teeth, irreversible damage to the neighboring teeth, breaking of the implant, prothesis, jaw or bone, esthetic problems, local or systemic allergic response. Some cases osseointegration could be failed.

It is the responsibility of the dentist to select patients who are expected to be satisfied psychologically, functionally and esthetically after recovery.

6. Warning

Implantation surgery and implantation prosthetics are complex dental procedures. It is especially recommended to master all necessary elements of proficiency. An improper surgical procedure may result in implant loss or bone degeneration. The implantation must be preceded by extensive diagnostic procedures. Improper design may result in unsuccessful implantation. Denti implants must only be loaded with devices included in the Denti Implant Systems and in accordance with guidelines and recommendations supplied therewith.

Products from different Denti Implant Systems are incompatible with each other, only abutments for implants from the same product system can be used. If you are using the DN and DNB implants, a dental check-up every six months is recommended.

Studying the instructions for use alone is not enough to start utilizing Denti Implant Systems. Denti implants may exclusively be loaded by oral surgeons, dentists and specialists trained to use Denti Systems. Denti System Ltd. holds regular trainings, which allow for a deeper knowledge in this field. Guarantees shall be submitted only if the claimer is able to confirm the completion of a Denti System Ltd. course.

7. Sterilization and packing

Implants are released following gamma sterilization and are made for a single use only. Clinically infected implants must not be cleaned and resterilized. Denti System Ltd. does not take any liability for the utilization of Denti implants taken from open, resterilized Denti implant sets. Used or unsterile implants must not be loaded. The sterile packing protects Denti implants from outer contamination. Inner packing must only be opened directly before loading. Prior to loading, personnel should check that the packing was not damaged in a way that the implant sterility was compromised. When stored properly (at a dry place, protected from sunlight at room temperature and in its original packing) implants remain sterile till the date of expiry printed on the packing. Implants must not be used after the date of expiry!</

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